

Extraction

12/1/08

56

Name

Course

Section

Lab partner

chem
226







3A

+10

Objective: Use solvent extraction techniques to separate a mixture consisting of a carboxylic acid, a phenol, and a neutral compound.

GS
GS

Physical Properties:

Compound Name	Structure	CAS #	BP or MP lit °C	RI	Hazards
Acetanilide		103-84-4	BP: 304 MP: 111-115	N/A	Irritating, harmful if swallowed
P-tert-butyl phenol		98-54-4	BP: 236-238 MP: 96-100	N/A	Irritant, harmful to the environment
Hydrochloric Acid	H-Cl	7647-01-0	BP: 57 MP: -35	N/A	Irritant, corrosive
Methyl-tert-butyl ether		1634-04-4	BP: 55-56 MP: -109	1.3675- 1.3695	Highly flammable, Irritant
Sodium Bicarbonate		144-55-8	BP: N/A MP: 270	N/A	Minor irritant
Sodium Hydroxide	Na ⁺ HO ⁻	1310-73-2	BP: 1390 MP: 318	N/A	Corrosive, causes burns, exothermic in contact w/H ₂ O
Sodium Sulfate		7757-82-6	BP: 1700 MP: 884	N/A	Irritant
P-toluic acid		99-94-5	BP: 274-275 MP: 179-182	N/A	Irritant, harmful to eyes, skin, resp. system

References:

- 1.) 30 Nov 2008 <<https://itweb.mst.edu/auth-cgwrap/msd.shtml/search.pl>>
- 2.) 30 Nov 2008 <<http://www.chemexper.com/>>

Procedures:

- 1.) A centrifuge tube was obtained.
- 2.) 5 mL of ~~the~~^{the} sample was added.
- 3.) 2 mL of 0.5M aqueous NaHCO_3 was added to the ether solution in the centrifuge tube.
- 4.) The two layers in the centrifuge tube were gently and thoroughly mixed.
- 5.) The plastic cap was placed on the centrifuge tube and shaken gently.
- 6.) The cap was removed to allow any CO_2 gas to escape.
- 7.) This process was repeated several times with the shaking intensity gradually increasing.
- 8.) The centrifuge tube was supported in a beaker and the layers were allowed to separate.
- 9.) The identity of the layers was confirmed by using a Pasteur pipet to introduce one or two drops of water just below the surface of the top layer. The aqueous layer was the bottom layer.
- 10.) A Pasteur pipet was used to remove the aqueous layer and transfer it to a 50 mL beaker.
- 11.) A second 2 mL NaHCO_3 portion was added to the tube containing the ether mixture.
- 12.) The tube was shaken vigorously.
- 13.) The aqueous layer was removed and combined in the beaker with the first extract.
- 14.) A third 2 mL NaHCO_3 portion was used to repeat the process.
- 15.) 1 mL of distilled water was added to the centrifuge tube and mixed.
- 16.) The aqueous layer was removed and combined with the three NaHCO_3 solution extracts in the 50 mL beaker.
- 17.) The aqueous solution was saved for later.
- 18.) 2 mL of 0.5M NaOH was added to the ether solution remaining in the centrifuge tube. The tube was shaken vigorously.
- 19.) A Pasteur pipet was used to remove the aqueous layer (bottom) and transfer it to a clean 50 mL beaker.
- 20.) A second 2 mL NaOH portion was used to repeat the

Extraction

12/1/08

58

Course

Section

Lab partner

Chem 226 3A

extraction of the ether layer.

- 21.) The second NaOH layer was removed and combined with the first.
- 22.) This was repeated using a third 2mL NaOH portion.
- 23.) 1 mL of water was added to the ether remaining in the centrifuge tube and mixed.
- 24.) The aqueous layer was removed and combined with the three NaOH extracts.
- 25.) The NaOH extracts in the 50mL beaker was saved for later. The remaining ether layer in the centrifuge tube for later.
- 26.) The 50mL beaker containing the NaHCO_3 extracts was selected.
- 27.) 3 mL HCl was added dropwise to the NaHCO_3 solution.
- 28.) The 3M HCl was added dropwise with stirring until no more solid was produced and the solution tested acidic.
- 29.) The crystals from the solution were separated from the solution using vacuum filtration with a Hirsch funnel.
- 30.) The 50 mL beaker containing the NaOH extracts was selected.
- 31.) The NaOH solution was heated to 60°C on a hot plate.
- 32.) The beaker was then allowed to cool.
- 33.) 3M HCl was added dropwise to the cooled solution until it tested acidic.
- 34.) To facilitate crystallization the mixture was cooled in an ice bath and the sides and bottom of the beaker were ~~scratched~~ ^{scratched} with a microspatula.
- 35.) The crystals were separated from the solution by vacuum filtration using a Hirsch funnel.
- 36.) The crystals were allowed to air dry.
- 37.) The centrifuge tube containing the ether layer was selected.
- 38.) .5g anhydrous Na_2SO_4 was added to the centrifuge tube.
- 39.) The tube was capped, shaken and allowed to stand for 5min.
- 40.) A 50mL beaker was weighed 25.578g.
- 41.) The dried ether was decanted into the 50mL beaker, and the Na_2SO_4 was left in the centrifuge tube.
- 42.) The ether was evaporated on a hot plate at 50°C with air blowing gently over the solution.

Extraction

12/1/08

Course


Section

Lab partner

chem226

3A

- 43.) The beaker was cooled in an ice bath. The bottom and sides of the beaker were scratched with a microspatula.
- 44.) The crystals were allowed to dry.
- 45.) The mass of the 1st extract was .14 grams.
- 46.) The mass of the 2nd extract was .264 grams.
- 47.) The mass of the 3rd extract was .065 grams.
- 48.) All materials were cleaned and put away or properly disposed of.



Vacuum Filtration Apparatus

Observations:

- The first crystals obtained were white
- The second set of crystals were whitish brown
- The beaker needed to be scratched on the bottom and the sides to induce crystallization

extract number	yield in grams	Melting Point	Percent Yield	Melting point % error
extract 1	.14	170-172 °C	2%	5.5%
extract 2	.264	95-96 °C	.66%	4%
extract 3	.065	110-112 °C	4%	2.6%

Calculations: Percent Yield

$$\text{extract 1: p-toluic acid } 6.5\text{g} \times \frac{1\text{ mole}}{136.15\text{g/mol}} = .048\text{ mol}$$

$$.048\text{ mol} \times \frac{136.15}{1\text{ mole}} = 6.5352\text{ grams}$$

$$\frac{.14}{6.5352} \times 100\% = 2\%$$

$$\text{extract 2: p-tert butylphenol } 6.5\text{g} \times \frac{1\text{ mole}}{150.22} = .043\text{ mol}$$

$$.043\text{ mol} \times \frac{150.22}{1\text{ mol}} = 6.46\text{ grams}$$

$$\frac{.043}{6.46} \times 100\% = .66\%$$

$$\text{extract 3: Acetaminilide } 4\text{g} \times \frac{1\text{ mole}}{135.16} = .030\text{ mol}$$

$$.030\text{ mol} \times \frac{135.16}{1\text{ mol}} = 4\%$$

Melting Point Percent error

$$\text{extract 1 p-toluic acid } \frac{|171 - 181|}{181} \times 100\% = 5.5\%$$

$$\text{extract 2 p-tert butylphenol } \frac{|95.5 - 99.5|}{99.5} \times 100\% = 4\%$$

$$\text{extract 3 Acetaminilide } \frac{|111 - 114|}{114} \times 100\% = 2.6\%$$